



Tube-Fix Securement of endotracheal tubes



Tube-Fix

Securement of endotracheal tubes

Securement device for endotracheal tubes

Tube-Fix is especially designed to secure the position of oral endotracheal tubes for artificial ventilation with minimum risk of creating pressure points.

Two versions

Ref. No. 650M is designed for securement of oral endotracheal tubes together with guedel airways, and Ref. No. 655M for securement of oral endotracheal tubes alone.

Both Ref. Nos. 650M and 655M allow easy access to the mouth and pharynx for suction and oral care.

Two different materials

Tube-Fix is available in two types of materials. Both materials are manufactured with a skin-friendly medical grade adhesive, suitable for both the patients face and the oral endotracheal tube.

REF.	Description	Inner Box Quantity	Outer Box Quantity
650M-F	Endotracheal Tube & Guedal Airway Endotracheal Tube	100 pcs	2000 pcs
650M-M	& Guedal Airway	100 pcs	2000 pcs
655M-F 655M-M	Endotracheal Tube Endotracheal Tube	100 pcs 100 pcs	2000 pcs 2000 pcs

All Unomedical products are CE-marked according to the Directive 93/42/EEC concerning medical devices.

• Foam

For fixation in operating or emergency rooms where it may be exposed to liquids.

Microporous non-woven material
 A peach coloured material which allows the
 skin to breathe. The microporous material is
 seen in many intensive care departments,
 and is suitable for short term use.



REF 655M-F (Foam)

REF 655M-M (Non-wowen)

Material specifications

Foam: Adhesive on foam: Microporous material: Adhesive on micropore: Extra adhesive layer on PE: Backing paper: Pouch:

Shelf life:

In situ:

White Polyethylene foam Medical grade acrylic adhesive Non-woven polypropylene Medical grade acrylic adhesive Coated on both sides with medical grade acrylic Kraft paper siliconized on one side Polyethylene 2 years foam version

3 years micropore version 24 hours For further information and availability please contact:

Unomedical a/s Birkerød Kongevej 2 DK-3460 Birkerød Denmark

Phone: +45 48 16 70 00 Telefax: +45 48 16 70 45 www.unomedical.com